

Clinical trials of tPA+heparin

TrialResults-center www.trialresultscenter.org

1 venous thrombosis

Trial	Treatments	Patients	Trials design and methods
tPA+heparin vs no fibrinolysis			
Schweizer tPA , 1998 n=NA follow-up:	tPA 20 mg IV into pedal vein over 4 hours each day for 7 days. Heparin IV given concomitantly, with adjustment versus heparin IV, adjusted for 7 days	patients with venographically confirmed DVT of leg duration <7 days.	Parallel groups single blind Germany

More details and results :

- fibrinolysis for venous thrombosis in all type of patients at <http://www.trialresultscenter.org/go-Q100>

References

Schweizer tPA, 1998:

Schweizer J, Elix H, Altmann E, Hellner G, Forkmann L Comparative results of thrombolysis treatment with rt-PA and urokinase: a pilot study. *Vasa* 1998;27:167-71 [9747153]

Entry terms: heparin, Heparin, Unfractionated Heparin, Heparinic Acid, Liquaemin, Sodium Heparin, Heparin Sodium, alpha-Heparin, alpha Heparin,